

BMW Motorrad

Maintenance schedule

0477 - R 850 RT



00 00 119 BMW Running-in Check (see Maintenance Schedule) (once only at 1,000 km)

Customer _____ Registration No. _____ Odometer reading _____

Order No. _____ Date _____ Mechanic's + Manager's signature _____

- | | |
|---|--------------------------|
| Reading fault-code memories with the BMW Motorrad diagnostic system | <input type="checkbox"/> |
| Performing bleed test with the BMW Motorrad diagnostic system | <input type="checkbox"/> |
| Changing engine oil with engine at operating temperature, replacing oil filter | <input type="checkbox"/> |
| Changing oil in final drive, warmed to normal operating temperature | <input type="checkbox"/> |
| Checking swinging-arm pivot mount | <input type="checkbox"/> |
| Checking operation of side-stand switch | <input type="checkbox"/> |
| Greasing side stand pivot | <input type="checkbox"/> |
| Tightening cylinder head | <input type="checkbox"/> |
| Adjusting valve clearance | <input type="checkbox"/> |
| Checking freedom of movement of throttle cable and checking for kinks and chafing | <input type="checkbox"/> |
| Checking/adjusting synchronisation | <input type="checkbox"/> |
| Checking tyre pressures and tread depth | <input type="checkbox"/> |
| Checking lights and signalling equipment | <input type="checkbox"/> |
| Test ride as final inspection and function check | <input type="checkbox"/> |
| Confirming BMW Service in on-board documentation | <input type="checkbox"/> |

BMW Motorrad, UX-VA-4
18.05.06

BMW recommends